


Elastomer Half Sine Performance Kit SD16 & SD24 Series and Auto Shock-II





Pulse #	Elastomer Combination				
1	1 in. Red Open				
2	1 In. Blue Closed, .5 in. Red Open				
3	1 in. Red Open, .5 in. Red Open				
4	1 in. Blue Closed, 1 in. Blue Closed, 1 in. Red Open				
5	2 in. Blue Closed, 1 in. Blue Open				
6	1 in. Blue Closed, 1 in. Blue Closed, 2 in. Blue Closed, 1 in. Blue Open				
7	1 in. Blue Closed, 1 in. Blue Closed, 2 in. Blue Closed, 1 in. Red Open				
8	1 in. Blue Closed, 1 in. Blue Open				




The Above Curves are nominal and contingent upon

- Payload
- Machine Operational Condition
- Elastomer Wear
- Temperature






Pulse #	Elastomer Combination				
9	1 in. Blue Closed, 1 in. Blue Open				
10	1 In. Blue Open				
11	2 in. Blue Open				
12	1 in. Green Closed, 1 in. Blue Open				
13	2 in. Blue Open, 1in. Blue Open				
14	2 in. Green Closed, 2 in. Blue Open				
15	2 in. Green Closed, 1 in. Green Closed, 2 in. Blue Open				
16	2 in. Blue Closed, 1 in. Blue Closed, 2 in. Green Closed, 1 in. Green Closed, 2 in. Blue Open				



The Above Curves are nominal and contingent upon

- Payload
- Machine Operational Condition
- Elastomer Wear
- Temperature





Elastomer Kit- L.AB. Item # 604044 Includes the following elastomers modules.

L.A.B. Item Number	Qty	Type	Color	Size
400544	1	Open	Red	1 Inch
400545	1	Open	Blue	1 Inch
400546	1	Open	Blue	2 Inch
400547	1	Open	Red	1/2 Inch
400548	2	Closed	Blue	1 Inch
400549	1	Closed	Green	1 Inch
400550	1	Closed	Blue	2 Inch
400551	1	Closed	Green	2 Inch

Note: Size referes to the appoximate thickness of the elastomer material not the assmbled Item number.



YOUR ONLY SOURCE FOR CUSTOMER-FOCUSED SHOCK AND VIBRATION TEST EQUIPMENT

L.A.B. Equipment, Inc.